

# Unistat Grande Fleur

**Unistat Grande Fleur controls the process temperature in non-insulated 1l glass jacketed reactor from Buchi**



### Requirement

This case study demonstrates the ability of the Unistat Grande Fleur to control the process temperature in non-insulated 1l Buchi glass jacketed reactor.

### Method

The Unistat Grande Fleur was connected to a 1l Buchi non-insulated glass jacketed reactor via 2 x 1-meter metal insulated tubes. The HTF used was Huber M40.165/220.10 and the process mass simulated with 0,7l of Huber M40.165/220.10 silicon oil.

Under "Process Control" using a Pt100 (Teflon covered), located in the process mass, different set-points were entered and the performance of the Unistat Grande Fleur was recorded using Huber's software.

The agitator speed was set to 250rpm.

### Setup details


- Temperature range: -40°C...+200°C
- Heating power: 1.8 kW
- Cooling power: 0.6 kW @ +100°C
- 0.6 kW @ 0°C
- 0.35 kW @ -20°C
- 0.2 kW @ -30°C
- Hoses: 2 x M24x1m Metal Insulated
- HTF: M40.165/220.10
- Reactor: Buchi 1l
- Reactor content: 7l M40.165/220.10
- Control: process
- Stirrer speed: 250 rpm
- Amb. temperature: +25°C

## Results

### 1. Temperature Control: from +20°C to +100°C

This test demonstrates the speed and accuracy that the Unistat Grande Fleur controls the process temperature from +20°C to +100°C and back to +20°C.

Start (°C)	End (°C)	Approximate time (min)	Average Ramp Rate (K/Min)
+20	+100	20	4.00
+100	+20	28	2.86



## 2. Temperature Control: -20°C to +100°C to -20°C

This test demonstrates the speed and accuracy that the Unistat Grande Fleur controls the process temperature from -20°C to +100°C and cooling down to -20°C.

Start (°C)	End (°C)	Approximate time (min)	Average Ramp Rate (K/Min)
-20	+100	29	4.13
+100	-20	54	2.22

